













26/08/2024

## TECHNICAL DATASHEET

Code	BIO-3000-SB
<b>Brief Description</b>	UN 4GU/CLASS 6.2 (UN2814/UN2900) complete packaging
Print	Neutral
Provided	Assembled with 95kPa tested secondary container and various products
Approval	4GU/CLASS.6.2/** B/SERPAC.403-090064
In compliance with	P620 ADR, P620 RID, P620 IMDG, P.I. 602 IATA
Capacity	0,95 L
Certification	95kPa from -40 °C to +55 °C
Information	Compliant with the requests of the Ministry of Health for the transport of Coronavirus swabs (Circular no. 3 of 8 May 2003).
Country of origin	Italy
Dimensions inner	ø100 x 149 h mm
Dimensions outer	130 x 130 x 180 h mm
Package	of 16 pieces



For shipping:

<b>*</b>	ADR	•
✓	ADR	•









Approved for	Content
Primary container of all types (test-tubes, vials, flacon, appendorf)	Outer fibreboard packaging UN certified 4GU/Class 6.2, 95kPa tested from +40°C to +55°C secondary rigid plastic container, cardboard insertion, nr. 1 absorbent pouch with 6 slots, 1 bubble-wrap pouch with 6 slots, label LAB-T621-PC (Infectious Substance), instruction sheet, security seal. COVID-19 packaging suitable for shipping Coronavirus infected samples.

## Description

## **Packaging for Infectious Substances**

\*Supplied complete with: secondary container, cardboard insertion, 1 absorbent test-tube holder with 6 slots, 1 bubble-wrap pouch with 6 slots, LAB-T621-PC label (Infectious Substance), instruction sheet, security seals and outer cardboard packaging. Always available in stock packaging certified for the transport of Infectious Substance Category A. Infectious Substance Category A are

- UN 2814 INFECTIOUS SUBSTANCE AFFECTING HUMANS
- UN 2900 INFECTIOUS SUBSTANCE AFFECTING ANIMALS. The Infectious substances may, when exposure occurs, causing permanent disability or life-threatening or fatal disease to people or animals, which had been in good health. For shipping Infectious Substance Category A (UN2814 – UN2900) you have to use a packaging certified for this substance. Our

packagings are certified UN 4GU/Class 6.2 and in the rigid plastic secondary container can be inserted primary containers of any type (test tubes, blood vials, falcon, appendorf, swabs and sample containers).

This UN certified packaging can be used for testing Coronavirus swabs.